


Hygiene & Disinfection



Dedicated Healthcare Solutions



Infection Control

About Uhs

Established in 1978, Uhs has led the way in offering cost-effective solutions to healthcare providers for over 25 years.

Essentially a Sales & Marketing operation, Uhs has a vast portfolio of products that originate from:

- Its own manufacturing facility
- OEM products from UK and overseas manufacturers
- Exclusive partnerships with other major health care manufacturers: Kimberly-Clark, Sterisol UK, PDI, Safeguard, Exel, Antec International & Yamasu.

As one of the major suppliers to NHS Logistics, we utilise the services of:

- Southern Syringe Services Ltd, the largest distributor of disposable medical and surgical products to UK hospitals
- 3S Health Care, providing the Uhs product portfolio to all other areas of the healthcare market
- To order all Universal Products, please refer to details on back cover.

The Hygiene and Disinfection Selection

The Uhs Hygiene & Disinfection range of products have definite applications for infection control in healthcare in 3 distinct areas:

Hand Hygiene and disinfection / decontamination of reusable medical devices such as flexible endoscopes and ultrasound probes and the cleaning and disinfection of surfaces in the patient environment, clinical areas and community based healthcare settings. This comprehensive range of products has been specifically picked for their quality, safety, ease of use and cost effectiveness.

Contents

Product Selection

Sterisol Skin Care System	3
The Reinol Range	6
Antec/Dupont Sterilisation	7
PDI Wipes Range	10

A Guide to Using Some of Our Products

Sterisol Hand Hygiene System	
- Using Disinfectant	17
Sterisol Hand Hygiene System	
- Using Soap	18
RelyOn™ Virkon Tablets / Spray Bottles	
Preparation Instructions	19
RelyOn™ PeraSafe Powder	
Preparation Instructions	19

All Products & Code Information

All code information can be found on pages	20-23
--	-------

For further products & code information see page 20

STERISOL SKIN CARE SYSTEM

Sterisol are the leading Scandinavian manufacturer of modern skin care products, Universal offer a system to wash, disinfect and condition the skin.

The Sterisol System contains no preservatives or artificial colourings, thus reducing skin allergies and irritation. Presented in an airtight cartridge with a unique valve that ensures bacteria-proof dispensing and correct product dosage.

Sterisol also offer a range of specially designed wall brackets and holders to meet the individual requirements of all healthcare settings.

Hand Soap and Cleansers

Hand hygiene is widely acknowledged to be the single most important activity for reducing the spread of disease. Hands that are visibly dirty or contaminated must be washed with Soap or Hand Cleanser and then dried thoroughly.

Hand Cleanser - Cleans resistant soiling on the hands without damaging the skin's barrier function.

- Do not contain solvent or scrubbing particles
- Free from preservatives.

Hand Soap - A soap that is kind to the skin. Excellent for those who wash many times a day. Also available:

- Ultra Mild
- Unscented
- All products are free from preservatives.



1 Uhs Ref: UNS7810 ■ Liquid Soap - Unscented 700ml

■ Pack QTY: 12

2 Uhs Ref: UNS7481 ■ Hand Cleanser 700ml

■ Pack QTY: 12

For further products & code information see page 20

Hand Disinfectants

Proper hand hygiene and regular hand disinfection are the best ways of preventing the spread of infection. Sterisol hand disinfectants destroy the potentially harmful transient bacteria, they also contain glycerine which counteracts dehydration and keeps the skin in good condition. Sterisol hand disinfectant is available in 3 dispensing formats:

- 700ml wall dispensing units
- 700ml surface mounted dispensing globe
- 120ml pocket size bottle.

Available in Isopropanol and Ethanol formulations.



1 Uhs Ref: UNS7921 ■ Hand Disinfectant - Iso/Nproponol 700ml

■ Pack QTY: 12

2 Uhs Ref: UNS7006 ■ Hand Disinfectant - Globe (Ethanol) 700ml

■ Pack QTY: 6

Hand Creams and Lotions

Regular use of a skin cream or lotion will help prevent hands from becoming dry, dehydrated and chapped. When the skin surface becomes damaged or broken it becomes more susceptible to infection and irritation. Both Sterisol products will help moisturise the skin, the creams are richer than lotions and would be more suitable for hands that are frequently washed over a short period of time. Sterisol Creams and Lotions:

- Do not contain preservatives
- Free from artificial colourings
- Available in perfume free versions.



1 Uhs Ref: UNS7359 ■ Skin Lotion 700ml

■ Pack QTY: 12

2 Uhs Ref: UNS7343 ■ Skin Cream 700ml

■ Pack QTY: 12

For further products & code information see page 20

Industrial and Recreational

Shower Soap, Shampoo and Liquid Soap/ Shampoo Ultra Mild are ideal skin care products for the hospital shower and bath facilities. The products are conveniently packaged, suitable for both body and hair and therefore very easy to use when washing patients.

- Available unscented and Ultra mild versions
- Free from preservatives.



1 Uhs Ref: UNS4808 ■ Hair & Body 2 in 1 Soap 5 Litre

■ Pack QTY: 2

2 Uhs Ref: UNS4531 ■ Dispenser 5 Litre

■ Pack QTY: 4

Dispensers

Sterisol dispensers are designed for quick and easy installation. By using the same back plate the complete range of sterisol dispensers become interchangeable without having to re-drill unsightly holes in the wall. The back plate can be easily and securely fixed to any wall using self adhesive pads or screws provided. The dispensers also benefit from:

- **Easy to clean design** – Dispensers can be unclipped from the back plate and cleaned effectively
- **No maintenance** – Strong design moulded from tough plastic to withstand heavy use
- **Minimising cross contamination risks** – Each product bag contains its own bacterial tight dispensing valve. Non drip, non clog dispensing.



1 Uhs Ref: UNS4580 ■ Washroom Dispenser 700ml

■ Pack QTY: 1

2 Uhs Ref: UNS4581 ■ White Dispenser 700ml

■ Pack QTY: 6

For further products & code information see page 20

THE REINOL RANGE

Universal provide a full range of Reinol Hand Cleaners, Skin Care and Barrier Creams.

The Reinol product are an environmentally friendly and economical product range that have established new standards in occupational skin care and health and safety in the workplace.

REINOL is a NON SOLVENT Hand Cleaner formulation that lifts and absorbs the dirt into the Hand Cleaner, that is then easily washed away.

REINOL Hand Cleaner cleans the dirtiest of hands with only one teaspoon of cleaner and simultaneously nourishes and conditions the skin. Economical and Skin Friendly Reinol is entirely Solvent Free, made from 16 natural ingredients and essential oils coupled with Reinol's own unique quartz grime exfoliators. Can be used for medium to heavy duty applications.

1



2



3



4



1 Uhs Ref: UNR3717 ■ Reinol Dermasoft 100ml Tube

■ Pack QTY: 12

2 Uhs Ref: UNR5161 ■ Reinol Paste Dispenser

■ Pack QTY: 1

3 Uhs Ref: UNR5031 ■ Reinol S 10 Litre Bucket

■ Pack QTY: 1

4 Uhs Ref: UNR5144 ■ Reinol Economic 1 Litre

■ Pack QTY: 10

For further products & code information see page 20

ANTEC/DUPONT STERILISATION

DuPont Relyon™ products have applications for infection control in hospitals in 2 distinct areas:

- Decontamination of reusable medical devices such as flexible endoscopes and ultrasound probes

- Cleaning with disinfection of surfaces in the patient environment, clinical areas, equipment and decontamination of body fluid spills.



Perasafe Instrument Sterilant Powder

RelyOn™ PeraSafe™ is a safer form of peracetic acid which is rapidly effective for decontamination of medical devices having proven efficacy against Mycobacterium tuberculosis, H.I.V. and other blood-borne viruses as well as vegetative and spore forming bacteria. RelyOn™ PeraSafe is a rapid, medical device compatible, user and environmentally friendly chemical sterilant.

Supplied as a safe, inactive powder, RelyOn™ PeraSafe is activated by simply dissolving in tap water. Providing rapid decontamination of medical devices, PeraSafe™ gives confidence to users and patients alike. An aldehyde-free formulation gives efficacy without harm to equipment, health risks to staff and avoids the need for special fume extraction or waste disposal systems.

RelyOn™ PeraSafe powder is available in a range of packs for varying user requirements. For example: multiples of 81 gm/5 Litre for endoscope reprocessors, 16.2 gm/1 Litre sachets for small volume, manual reprocessing, for example, lenses in eye clinics, dental instruments and naso-endoscopes and instruments in Out Patients.

- | | | |
|---|---|-----------------|
| 1 | Uhs Ref: UND12554795 ■ Persafe Sterilant 810g 50 Litre | ■ Pack QTY: 1 |
| 2 | Uhs Ref: UND12554755 ■ Persafe Sterilant 81g 5 Litre | ■ Pack QTY: 24 |
| 3 | Uhs Ref: UND12554720 ■ Perasafe Sterilant Powder 16.2g 1 litre sachet | ■ Pack QTY: 100 |

For further products & code information see page 20

Perasafe Instrument Sterilant Liquid

Based on the same chemistry as the original powder RelyOn™ PeraSafe, the inactive and safe to handle components are in a binary liquid consisting of a Base and Activator.

Equal volumes of Base and Activator are diluted and mixed with water to form the active sterilant solution. The solution is at a neutral pH, which is optimal for compatibility and efficacy across the microbiological spectrum.

RelyOn™ PeraSafe Liquid is available in 2 x 5 Litre packs. The Multiple Use packs are supplied with dispense pumps to facilitate measuring. Solutions from 1 Litre upwards can be easily prepared. Packs generate 165 Litres of active solution in total.

RelyOn™ PeraSafe Liquid active solutions have been validated for efficacy against spores, mycobacteria, viruses, bacteria and fungi by recognised testing laboratories such as H.I.R.L. at City Hospital Birmingham.

Materials compatibility is assured through tests and years of experience with the original RelyOn™ PeraSafe.



1 Uhs Ref: UNA0673 ■ Persafe Multiple Use Liquid 2 x 5 Litre

■ Pack QTY: 1



Virkon Disinfectant

RelyOn™ Virkon is a reliable broad spectrum high level disinfectant which combines effectiveness, safety and ease of use with environmental compatibility.

RelyOn™ Virkon is extremely safe and easy to use. Ideal for maintaining a hygienic environment to reduce the transmission of pathogenic micro-organisms through contact with surfaces and equipment.

RelyOn™ Virkon® is especially effective for decontamination and prevention of spread of norovirus infections such as "Winter Vomiting Disease" and is effective for routine disinfection of environments where Methicillin Resistant Staphylococcus aureus gram positive bacteria can cause infections.

1 UND12555734 ■ Virkon Disinfectant 500g Shaker

■ Pack QTY: 6

2 UND12557258 ■ Virkon Tablets 5g 500ml

■ Pack QTY: 10

For further products & code information see page 20

Perasafe and Virkon Extras

This range of accessory products have been designed to facilitate the measuring, mixing and application of the Virkon and Perasafe range of products.




1




2



3



4



1	Uhs Ref: UND12564931 ■ Spray Bottle 500 ml	■ Pack QTY: 1
2	Uhs Ref: UND12564869 ■ Perasafe Test Strips 60 Per Pot	■ Pack QTY: 1
3	Uhs Ref: D12564953 ■ Plastic Mixing Jug 5 Litre	■ Pack QTY: 1
4	Uhs Ref: D12564961 ■ Stirring Rod 245mm Bag of 10	■ Pack QTY: 10

For further products & code information see page 20

PDI-WIPES RANGE

PDI is a leader in developing, testing and manufacturing pre-moistened wipes for infection control and prevention, education and patient care. Universal and PDI have created an infection control "umbrella" serving the healthcare market through hospitals, Doctors and dental offices, emergency medical services, long-term and home health care, blood banks and laboratories UK wide.

Personal Hygiene Wipes Hygēa Hands & Face Wipes Canister

- Available in 125 and 200 Wipe Canisters
- Large wipes offering high strength
- Hypoallergenic ensuring it's gentle to sensitive skin
- New slim-line 200 Wipe Canister is easier to handle, and has a Deep Well Lid providing better dispensing
- **Wipe size (cm):** 125 Canister = 18.5 x 13, 200 Canister = 20 x 22
- **Wipe weight (gsm):** 125 & 200 Canisters = 23
- **Units per case:** 125 & 200 Canisters = 6.



1 Uhs Ref: UNXP00161 ■ Hygēa Hands & Face Personal Wash Cloths 200 Wipes per Canister ■ Pack QTY: 6 Canisters

Personal Hygiene Wipes Hygēa Hands & Face Wipes Flow-wrap

- Available in 24 and 80 Wipe Packs and 80 Wipe Tub
- Wipes being pulp based are biodegradable
- Maceratable wipes ensure convenient and safe disposal
- Label offers space for patient/ward details
- High quality 50 gsm, offers both patient comfort & strength
- Lid on Tub, allows for one wipe to be dispensed at a time, without compromising the other wipes.
- **Wipe size (cm):** 17.5 x 20.3, **Wipe weight (gsm):** 50
- **Units per case:** 24 Pack = 24, 80 Pack = 12, 80 Tub = 6.



1 Uhs Ref: UNXP00166 ■ Hygēa Hands & Face Personal Wash Cloths 80 Wipes per Tub ■ Pack QTY: 6 Tubs

2 Uhs Ref: UNXP00164 ■ Hygēa Hands & Face Personal Wash Cloths - Patient Packs 24 Wipes per Pack ■ Pack QTY: 24

For further products & code information see page 20

Hard Surface Disinfectant Wipes Sani Cloth 70 Alcohol Wipes

- Available in 125 and 200 Wipe Canisters and as 100 count sterile disinfectant sachets
- Products exceed European Bactericidal standard - EN1276 with a fast initial kill of less than a minute
- Contains 70% alcohol

Canisters

- Canisters are pilfer sealed ensuring no evaporation and a quality product
- Deep well dispenser lid on 200 Wipes Canister allows for more controlled dispensing and less wastage

Sachets

- Sachet Wipes are sterile & gamma irradiated
- Sachets packed small for convenience, but the area of the wipe remains large.



1 Uhs Ref: UNXP00159 ■ Sani 70 Alcohol Wipes 200 Wipes per Canister

■ Pack QTY: 6 Canisters

2 Uhs Ref: UNWD24100 ■ Sani 70 Alcohol Wipes Individual Wrapped Sachets 100 Sachets per Box

■ Pack QTY: 1 Box

Hard Surface Disinfectant Wipes Super Sani Cloth Plus Wipes



- Available in a 200 Wipe Canister and as 100 count sterile disinfectant sachets
- Products exceed European Bactericidal standard - EN1276 with a fast initial kill of less than a minute
- Contains 70% alcohol and 1% quat
- **Wipe size (cm):** 19.6 x 12.7
- **Wipe weight (gsm):** 23
- **Units per case:** 10.

1 Uhs Ref: UNWD23800 ■ Super Sani Cloth Plus Wipes Individual Wrapped Sachets 100 Sachets per Box

■ Pack QTY: 1 Box

For further products & code information see page 20

Hard Surface Disinfectant Wipes Sani Cloth Active Wipes

- Alcohol Free wipes with a high spectrum of disinfection
- Canister is pilfer sealed ensuring no evaporation and a quality product
- Deep well lid provides greater protection of the wipe, and reduces wastage and contamination
- Meets European testing standards as Bactericidal – EN1040, EN1276, EN13697, Fungicidal – EN1275, EN1650, Viral – EN 14476, Polio, HepB, Herpes, Norwalk, Adenovirus and TB
- Product has also been shown to be effective against Bird Flu
- Product has been independently validated by BSI as being CE marked as a class 11a medical device, as required under rule 15 of the Medical Device Directive
- **Wipe size (cm):** 125 Canister = 13 x 22, 200 Canister = 20 x 22
- **Wipe weight (gsm):** 125 & 200 Canisters = 23
- **Units per case:** 125 & 200 Canisters = 6.



1 Uhs Ref: UNXP00188 ■ Sani Cloth Active Wipes 125 Wipes per Canister

■ Pack QTY: 6 Canisters

2 Uhs Ref: UNXP00190 ■ Sani Cloth Active Wipes 200 Wipes per Canister

■ Pack QTY: 6 Canisters

Hard Surface Disinfectant Wipes Sani Cloth Detergent Wipes



- Available in 225 count Detergent Wipe Bucket and 60 count flow-wrap – disinfectant free
- High gsm wipe offers better quality, for better cleaning results
- Handle on container is made of plastic offering hand protection
- Alcohol free
- Reduced foam
- For general damp dusting and routine cleaning of surfaces sensitive to alcohol
- Flow-wrap product is maceratable
- **Wipe size (cm):** 225 Bucket = 30 x 25, Flow Wrap 60 Pack = 30 x 20
- **Wipe weight (gsm):** 225 Bucket = 25, 60 Pack = 28
- **Units per case:** 225 Bucket = Single, 60 Pack = 9.

1 Uhs Ref: UNXP00175 ■ Sani Cloth Detergent Wipes 225 Wipes per Bucket

■ Pack QTY: 1 Bucket

2 Uhs Ref: UNDXP00186 ■ Sani-Detergent Flow-Wrap - Maceratable 60 Wipes per Pack

■ Pack QTY: 1 Pack

For further products & code
information see page 20

Disinfectant Swabs

Uni-Wipe Alcotip Swabs with 70% IPA

- Each swab individually sealed in it's own foil sachet
- 100 sachets per box.
- **Wipe size (cm):** 3 x 6.5
- **Wipe weight (gsm):** 45
- **Units per case:** 40.



1 Uhs Ref: UNG602 ■ Alcotip - 70% Alcohol Swab 100 per Box ■ Pack QTY: 20 Boxes

Oral Care



Lemon Glycerin Swabsticks 3's

Provides oral lubrication and stimulates saliva. Effective for patients with parched lips, on oxygen, after general anaesthesia and drug therapy. Also for patients prescribed fluid restricted diets. Lint free swabs are impregnated with Lemon Glycerin solution.

Saliswab Oral Care Swabstick 3's

Alleviates Xerostomia(dry mouth). May be used as a saliva substitute. Each Swab is saturated with a special formula that provides oral lubrication for patients on oxygen, general anaesthesia and drug therapy. Provides relief from the discomfort of dry mouth.

1 Uhs Ref: UNW16225 ■ PDI Lemon & Glycerin Swabstick Packed in 3s 25 x 3 per Box ■ Pack QTY: 10 Boxes
 2 Uhs Ref: UNW63925 ■ Saliswab Blackcurrant Swabstick Packed in 3s 25 x 3 per Box ■ Pack QTY: 10 Boxes

For further products & code information see page 20

Liquid & Ointments

Lubricating Jelly 2.7 gm

- Each packet terminally sterilized to prevent cross-contamination
- Makes routine procedures easier-catheter insertions and specula insertions
- Will not get tacky-maximizes patient comfort
- Water soluble-easily washed off. Contains Glycerin, Methyl and Propylparabens.



1 Uhs Ref: UNW00128 ■ PDI Lubricating Jelly 2.7gm 144 per Box

■ Pack QTY: 12 Boxes

Speciality Wipes Products

Electrode Skin Prep Pads

Cleanses skin in preparation for electrode placement. Pumice removes excess skin oil for better electrode-to-skin contact and longer interface-free contact. Contains 70% Isopropyl Alcohol and pumice.

Adhesive Tape Remover Pads

Effectively softens and removes tape and residue from intact skin. Reduces patient discomfort.

1 Uhs Ref: UNW59800 ■ Electrode Skin Prep Pads 100 Pads per Box

■ Pack QTY: 10 Boxes

2 Uhs Ref: UNW16400 ■ Adhesive Tape Remover Pads 100 Pads per Box

■ Pack QTY: 10 Boxes

For further products & code information see page 20

Speciality Wipes Products Continued

See Clear® Eye Glass Cleaning Wipes

Eyeglass cleaning wipe. Safely cleans eye protection gear and glasses. For glass and polycarbonate lenses. Anti-static / anti-fog formula. Dries quickly. Reduces eye strain. Does not scratch lenses. Packaged in convenient dispenser box.

Nail Polish Remover Pads

For the removal of nail polish when respiratory evaluation is necessary. Does not contain acetone. Contains Ethyl Acetate, Isopropyl Alcohol, Water, Propylene Glycol, Dicaprylate/Dicaprate, Butylene Glycol.



1 Uhs Ref: UNW19831 ■ See Clear Eyeglass Cleaning Wipes 120 Wipes per Box

■ Pack QTY: 12 Boxes

2 Uhs Ref: UNW71200 ■ Nail Polish Remover Pads 100 Pads per Box

■ Pack QTY: 10 Boxes

For further products & code information see page 20

Towelettes

Castile Soap Towelette

Ideal for midstream clean-catch urine specimen procedures. Gentle cleansing agent, pH balanced. Contains 2% Coconut oil-based soap solution.

Cleansing Towelette (5% alcohol)

Effective and gentle cleansing. Soothing and refreshing. Ideal for personal hygiene. Saturated with Benzalkonium Chloride 0.50%, Alcohol 5%.



1 Uhs Ref: UNW41900 ■ Castile Soap Towelette 100 Wipes per Box

■ Pack QTY: 10 Boxes

2 Uhs Ref: UNW31000 ■ Cleansing Towelettes - Individually Foiled Wrapped 100 Wipes per Box

■ Pack QTY: 10 Boxes

Wipes Wall Bracket



- The Wipes Wall bracket is designed for use with 200 count deep well canisters
- The Bracket is easy to install and maintain
- Will fit products UNXP00190, UNXP00161, UNXP00159.

1 Uhs Ref: UNXQ53521 ■ Wire Wall Bracket - Will fit UNXP00159, UNXP00141 & UNXP00161

■ Pack QTY: 1 Bracket

For further products & code information see page 20

A Guide to Using Some of Our Products


We are always striving to ensure that our products are used safely and correctly.

The following pages show the procedures for using some of the products detailed in the brochure, please read carefully and use as reference when needed. If you need further advice, contact your Uhs sales representative.

Products detailed: ■ Sterisol Hand Hygiene System - Disinfectant & Soap
■ Reylon™ Virkon Tablets/Spray Bottles ■ Reylon™ PeraSafe Powder

STERISOL HAND HYGIENE SYSTEM

USING DISINFECTANT




The technique shown is approved by:

ICNA - Infection Control Nurses Association*

*Ref: Gould.D. (1994) making sense of hand hygiene *Nursing Times*: 90:30, 63-64.

- most frequently missed
- frequently missed
- less frequently missed

 <p>1</p>	 <p>2</p>	 <p>3</p>	 <p>4</p>
Take a measure of disinfectant.	Work into hands, rub palm to palm.	Interlace fingers and rub.	Cup fingers and rub side to side.
 <p>5</p>	 <p>6</p>	 <p>7</p> <p>Make sure you rub in the solution until the hands are dry.</p>	
Rotate thumbs in palms of hands.	Rotate wrist in palms of hands.		

The above sequences are only a guide. Refer to your local hospital policies for further information.


STERISOL HAND HYGIENE SYSTEM

USING SOAP

The technique shown is approved by:

ICNA - Infection Control Nurses Association*

*Ref: Gould.D. (1994) making sense of hand hygiene *Nursing Times*: 90:30, 63-64.



- most frequently missed
- frequently missed
- less frequently missed

<p>1</p> <p>Wet hands and take a measure of soap.</p>	<p>2</p> <p>Work into hands, rub palm to palm.</p>	<p>3</p> <p>Interlace fingers and rub.</p>	<p>4</p> <p>Cup fingers and rub side to side.</p>
<p>5</p> <p>Rotate thumbs in palms of hands.</p>	<p>6</p> <p>Rotate wrist in palms of hands.</p>	<p>7</p> <p>Make sure you rinse all of the soap off your hands under running water and then dry your hands thoroughly with a paper towel.</p>	

Make sure you use hand cream after increased frequency in hand washing, before breaks and at the end of a shift, and remember... **Blue to clean, Red in between, Brown = Handcream.**



Blue = Soap



Red = Hand Disinfectant



Brown = Hand Cream

The system is colour coded, helping you to use the products in the correct order - use the colour guide above.

The above sequences are only a guide. Refer to your local hospital policies for further information.

RELYON™ VIRKON TABLETS / SPRAY BOTTLES PREPARATION INSTRUCTIONS

1 Fill container (bottle, beaker or spray bottle) with 500ml of tepid tap water.



2 Add 1 RelyOn™ Virkon tablet to the container and allow to dissolve.



3 Allow 5 minutes for tablets to dissolve and ensure solution is mixed thoroughly by stirring or agitation.



4 Solution is now ready to use.



RELYON™ PERASAFE POWDER PREPARATION INSTRUCTIONS



1 Fill a container with lukewarm tap water (30-35°C) to the required volume.



2 Add powder to water direct from container or required amount using measuring cap provided.



3 Stir until the powder dissolves to form a clear blue solution. (10-15 minutes).



4 Pour the solution into an immersion bath or endoscope reprocessor for use. The liquid will fade to a colourless solution in approximately one hour but can be used for up to 24 hours after activation. i.e. when the clear blue solution is achieved.

The above sequences are only a guide. Refer to your local hospital policies for further information.

Solution should be discarded if the test strip indicates it is below the Minimum Effective Concentration.

All Products & Code Information

Uhs supply a vast range of hygiene and disinfection equipment which is available in an array of sizes, quantities and colour references. It would be impossible to display all of these products in the catalogue; the product pages show part of the range.

The following pages detail the entire Uhs Hygiene and Disinfection range, together with code information, a small description and sizes, types and quantities available. If you wish to order products, note down the relevant codes and use the contact details on the back page, if you need further advice please contact your Uhs sales representative who will be pleased to help.

Code	Description		Qty
SECTION 2 - HYGIENE AND DISINFECTION PRODUCTS			
STERISOL SKIN CARE SYSTEM			
HAND SOAP AND CLEANSERS		Measure	
■ UNS7811	Liquid Soap	700ml	12
■ UNS7810	Liquid Soap - Unscented	700ml	12
■ UNS7815	Liquid Soap - Ultra Mild	700ml	12
■ UNS7816	Liquid Soap - Deluxe	700ml	12
■ UNS7842	Cleansing Cream	700ml	12
■ UNS7481	Hand Cleanser	700ml	12
■ UNS2000	Liquid Soap - Globe	700ml	6
■ UNS4486	Hand Cleanser - Ultra Clean	250ml	6
HAND DISINFECTANTS		Measure	
■ UNS7921	Hand Disinfectant - Iso/Nproponol	700ml	12
■ UNS7106	Hand Disinfectant - Ethanol	700ml	12
■ UNS7006	Hand Disinfectant - Globe (Ethanol)	700ml	6
■ UNS4110	Hand Disinfectant - Ethanol	120ml	24
■ UNS4120	Hand Disinfectant - ISO/Nproponol	120ml	24
HAND CREAMS AND LOTIONS		Measure	
■ UNS7343	Skin Cream	700ml	12
■ UNS7345	Skin Cream - Unscented	700ml	12
■ UNS7359	Skin Lotion	700ml	12
■ UNS7365	Skin Lotion - Unscented	700ml	12
INDUSTRIAL AND RECREATIONAL		Measure	
■ UNS7801	Shower Soap / Shampoo	700ml	12
■ UNS2001	Shower Soap - Globe	700ml	6
■ UNS4440	Hand Cleanser	5 Litre	2

For further product information
use the contact details below

Code	Description		Qty
■ UNS4484	Hand Cleanser Ultra	2.5 Litre	2
■ UNS4805	Shower Soap	5 Litre	2
■ UNS4808	Hair and Body - 2 in one Soap	5 Litre	2
■ UNS4814	Liquid Soap	5 Litre	2
■ UNS4817	Liquid Soap - Ultra Mild	5 Litre	2

DISPENSERS		Measure	
■ UNS4547	Dispenser - Plastic Arm & Transparent Cover	700ml	20
■ UNS4561	Dispenser - Transparent Cover	700ml	20
■ UNS4567	Dispenser - Metal Arm & Transparent Cover	700ml	20
■ UNS4569	Stainless Steel Dispenser	700ml	1
■ UNS8002	Dispenser - Globe	700ml	12
■ UNS8710	System 2000 Pole Adapter - Globe	700ml	12
■ UNS8706	Surface Stand for Globe	700ml	12
■ UNS4531	Dispenser	5 Litre	4
■ UNS4580	Washroom Dispenser	700ml	1
■ UNS4581	White Dispenser	700ml	6

REINOL RANGE			
■ UNR3717	Reinol Dermasoft 100ml tube		12
■ UNR5006	Reinol Original 1kg tub		24
■ UNR5007	Reinol Original 2kg tub		12
■ UNR5008	Reinol Original 3kg (2 litre) cartridge		6
■ UNR5011	Reinol Original 15kg (10 litre) bucket		1
■ UNR5024	Reinol S 500ml pot		30
■ UNR5028	Reinol S 2 litre cartridge		6
■ UNR5031	Reinol S 10 litre bucket		1
■ UNR5063	Reinol K 500ml pot		30
■ UNR5064	Reinol K 2 litre cartridge		6
■ UNR5065	Reinol K 3 litre cartridge		1
■ UNR5090	Reinol KMV 2.5 litre bottle		6
■ UNR5099	Reinolen 100ml tube		12
■ UNR5100	Reinol P Cream 100ml tube		12
■ UNR5116	Reinol Liquid HD 10 litre container		1
■ UNR5144	Reinol Economic 1 litre		10
■ UNR5566	Ja-Top 38 10 litre container		1
■ UNR5161	Reinol Paste Dispenser		1
■ UNR9998	Reinol Liquid Dispenser		1
■ UNR194259	Reinol S Dispenser		1
■ UNR194198	Reinol KMV Bracket		1
■ UNR194216	Reinol KMV Pump		1

Dispenser Spare Parts Available On Request

For further product information
use the contact details below

Code	Description		Qty
PERASAFE INSTRUMENT STERILANT		Weight (Measure)	
■ UND12554755	Perasafe Low Foam Sterilant	81g (5 Litre)	24
■ UND12554795	Perasafe Sterilant	810g (50 Litre)	1
■ UND12554720	Perasafe Sterilant	16.2g (1 litre sachet)	100
■ UND12554255	Perasafe Single Use Liquid	(2 x 5 Litre)	1
■ UND12554247	Perasafe Multiple Use Liquid	(2 x 5 Litre)	1
PERASAFE EXTRAS		Measure	
■ D12564953	Plastic Mixing Jug	5 Litre	1
■ UND12564869	Powder Perasafe Test Strips 60 Per Pot		1
■ D12564961	Stirring Rod 245mm Bag of 10		10
VIRKON DISINFECTANT		Weight (Measure)	
■ UND12555734	Virkon Disinfectant	500g Shaker	6
■ UND12555728	Virkon Disinfectant	5kg Drum	1
■ UND12555712	Virkon Disinfectant	50g Sachets	50
■ UND12557249	Virkon Tablets	5g (500ml)	50
■ UND12557258	Virkon Tablets	5g (500ml)	10
VIRKON EXTRAS		Measure	
■ UND12564931	Spray Bottle	500ml	1
■ D12564961	Stirring Rod 245mm Bag of 10		10
UNI-WIPE - WIPES RANGE			
ALCOHOL SWABS		Qty Per Box	
■ UNG602	Alcotip - 70% Alcohol Swab	Individually Foil Wrapped - 100 per box	1 box
PDI WIPES			
ORAL CARE		Qty Per Box	
■ UNW16225	PDI Lemon & Glycerin Swabstick Packed in 3s	25 x 3 per Box	10 Boxes
■ UNW63925	Saliswab Blackcurrant Swabstick Packed in 3s	25 x 3 per Box	10 Boxes
LIQUID & OINTMENTS		Qty Per Box	
■ UNW00128	PDI Lubricating Jelly 2.7gm	144 per Box	12 Boxes
■ UNW00229	PDI Lubricating Jelly 2.7gm	72 per Box	12 Boxes
■ UNW37370	PDI PVP Iodine Ointment 1.5gm	250 per Box	4 Boxes
TOWELETTES		Qty Per Box	
■ UNW41900	Castile Soap Towelette	100 Wipes per Box	10 Boxes
■ UNW31000	Cleansing Towlettes - Individually Foiled Wrapped	100 Wipes per Box	10 Boxes

For further product information
use the contact details below

Code	Description		Qty
PERSONAL HYGIENE WIPES		Qty Per Item	
■ UNW00295	Just Clean Skin Body Cloths - Individually Wrapped	255 Cloths per Box	1 Box
■ UNW12095	Hygea Wipes Premoistened Washcloth - Individually Wrapped	400 Cloths per Box	1 Box
■ UNXP00161	Hygea Hands and Face Personal Wash Cloths	200 Wipes per Canister	6 Canisters
■ UNXP00162	Hygea Hands and Face Personal Wash Cloths	125 Wipes per Canister	6 Canisters
■ UNXP00164	Hygea Hands and Face Personal Wash Cloths - Patient Packs	24 Wipes per Pack	24 Packs
■ UNXP00165	Hygea Hands and Face Personal Wash Cloths - Wipes Refill	80 Wipes per Pack	12 Packs
■ UNXP00166	Hygea Hands and Face Personal Wash Cloths	80 Wipes per Tub	6 Tubs
HARD SURFACE DISINFECTANT WIPES		Qty Per Item	
■ UNXP00142	Sani 70 Alcohol Wipes	125 Wipes per Canister	6 Canisters
■ UNXP00159	Sani 70 Alcohol Wipes	200 Wipes per Canister	6 Canisters
■ UNWD24100	Sani 70 Alcohol Wipes Individual Wrapped Sachets	100 Sachets per Box	1 Box
■ UNXP00188	Sani Cloth Active Wipes	125 Wiper per Canister	6 Canisters
■ UNXP00190	Sani Cloth Active Wipes	200 Wiper per Canister	6 Canisters
■ UNWD23800	Super Sani Cloth Plus Wipes Individual Wrapped Sachets	100 Sachets per Box	1 Box
DETERGENT WIPES		Qty Per Item	
■ UNXP00175	Sani Cloth Detergent Wipes	225 Wipes per Bucket	1 Bucket
■ UNDXP00186	Sani-Detergent Flow Wraps - Maceratable	60 Wipes per Pack	1 Pack
WALL BRACKET			
■ UNXQ53521	Wire Wall Bracket - Will fit UNXP00159, UNXP00141 and UNXP00161		1 Bracket
SPECIALITY WIPES PRODUCTS		Qty Per Box	
■ UNW16400	Adhesive Tape Remover Pads	100 Pads per Box	10 Boxes
■ UNW19831	See Clear Eyeglass Cleaning Wipes	120 Wipes per Box	12 Boxes
■ UNW59800	Electrode Skin Prep Pads	100 Pads per Box	10 Boxes
■ UNW71200	Nail Polish Remover Pads	100 Pads per Box	10 Boxes

Also available is the *Universal Equipment for Personal Protection Brochure*

This Product Selection includes:

Standard Face Masks, Fluid Resistant Face Masks, Speciality Face Masks, Protective Eye Wear, Moldex Respirators, Aprons, Head Wear Range and Gloves



**For further product information
use the contact details below**

A majority of Uhs products are available from the NHS Logistics Catalogue and Logistics Online.

Orders should be placed via our logistical partners. For all hospital orders contact your local depot:

Southern Syringe Services Ltd - London

George House, Unit 6 Delta Park Industrial Estate,
Millmarsh Lane, Enfield, London, EN3 7QJ
Telephone: 0870 873 4900 ■ Facsimile: 0870 873 4902
e-mail: sales@southern-syringe.co.uk

Southern Syringe Services Ltd - Manchester

Unit B Circle Court, Trafford Park,
Warren Bruce Road, Manchester, M17 1LB
Telephone: 0870 873 4918 ■ Facsimile: 0870 873 4919
e-mail: sales@manchester.southern-syringe.co.uk

Southern Syringe Services Ltd - Bristol

Enterprise House, Bradley Road,
Royal Portbury Dock, Bristol BS20 7XN
Telephone: 0870 873 4911 ■ Facsimile: 0870 873 4910
e-mail: sales@bristol.southern-syringe.co.uk

Southern Syringe Services Ltd - Leicester

New Frontier House, Interlink Park,
Bardon Hill, Coalville, Leicester LE67 1LE
Telephone: 0870 873 4916 ■ Facsimile: 0870 873 4917
e-mail: sales@leicester.southern-syringe.co.uk

Surgical Supplies Services Ltd - Scotland

St James House, 2 Tollpark Place, Unit 2 Wardpark
East Industrial Estate, Cumbernauld G68 0LN
Telephone: 0870 873 4913 ■ Facsimile: 0870 873 4912
e-mail: sales@cumbernauld.southern-syringe.co.uk

Macartney Surgical Ltd - Northern Ireland

Unit 3, Curran Business Park, Portland Road,
Larne, Co. Antrim BT40 1DH
Telephone: 0870 873 4914 ■ Facsimile: 0870 873 4915
e-mail: sales@larne.southern-syringe.co.uk

Non-Hospital Purchases:

3S Healthcare Ltd

George House, Unit 6 Delta Park Industrial Estate,
Millmarsh Lane, Enfield, London, EN3 7QJ
Telephone: 0870 873 4901 ■ Facsimile: 0870 873 4905



If you require any further product information contact our customer services team:

Universal Hospital Supplies

Unit 6, George House, Delta Park Industrial Estate,
Millmarsh Lane, Enfield, EN3 7QJ

■ Tel: 0845 082 0182 ■ Fax: 0845 082 0180
■ Email: sales@uhs.co.uk ■ Web: www.uhs.co.uk
■ Email: customerservices@uhs.co.uk

Dedicated Healthcare Solutions



Infection Control